

FMT

Medical Sensors, Cables and Accessories



ELECTROSURGERY CABLES & ACCESSORIES



Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

Generator Side Connector

Instrument Side Connector



Euro Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AA (3.0 m), ESU-BP/AAL (5.0 m)



Euro Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AAT (3.0 m), ESU-BP/AATL (5.0 m)



U.S. Forceps Connector (Angled)

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AB (3.0 m), ESU-BP/ABL (5.0 m)



U.S. Forceps Connector (Angled)

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/ABT (3.0 m), ESU-BP/ABTL (5.0 m)



U.S. Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş


REF: ESU-BP/AC (3.0 m), ESU-BP/ACL (5.0 m)



U.S. Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/ACT (3.0 m), ESU-BP/ACTL (5.0 m)



Aesculap Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AD (3.0 m), ESU-BP/ADL (5.0 m)



Aesculap Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/ADT (3.0 m), ESU-BP/ADTL (5.0 m)



R.Wolf MIS Instrument Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AE (3.0 m), ESU-BP/AEL (5.0 m)



R.Wolf MIS Instrument Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AET (3.0 m), ESU-BP/AETL (5.0 m)



K.Storz M.I.S. Instrument Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AF (3.0 m), ESU-BP/AFL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

Instrument Side Connector ▼



● **K.Storz M.I.S. Instrument Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AFT (3.0 m), ESU-BP/AFTL (5.0 m)



● **Ethicon/Everest Forceps Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş


REF: ESU-BP/AG (3.0 m), ESU-BP/AGL (5.0 m)



● **Ethicon/Everest Forceps Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AGT (3.0 m), ESU-BP/AGTL (5.0 m)



● **Ethicon Power Star Scissors Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AH (3.0 m), ESU-BP/AHL (5.0 m)



● **Aesculap Bipojet Scissors Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AI (3.0 m), ESU-BP/AIL (5.0 m)



● **Artery Sealer Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AJ (3.0 m), ESU-BP/AJL (5.0 m)



● **BiTech Scissors Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AK (3.0 m), ESU-BP/AKL (5.0 m)



● **To-Bite & Calvian Forceps & Coagulation Clamps Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AL (3.0 m), ESU-BP/ALL (5.0 m)



● **K.Storz Resectoscope Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AM(3.0 m), ESU-BP/AML (5.0 m)



● **Olympus HiQ+Hand Instrument Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/AN(3.0 m), ESU-BP/ANL (5.0 m)



● **Olympus HiQ+Hand Instrument Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş

REF: ESU-BP/ANT(3.0 m), ESU-BP/ANTL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

Generator Side Connector

Instrument Side Connector



Tontarra, Wolf Resectoscope Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf, Alan Combi Max

REF: ESU-BP/AO (3.0 m), ESU-BP/AOL (5.0 m)



Euro Forceps Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BA (3.0 m), ESU-BP/BAL (5.0 m)



U.S. Forceps Connector (Angled)

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BB (3.0 m), ESU-BP/BBL (5.0 m)



U.S. Forceps Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BC (3.0 m), ESU-BP/BCL (5.0 m)



Aesculap Forceps Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BD (3.0 m), ESU-BP/BDL (5.0 m)



R.Wolf MIS Instrument Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BE (3.0 m), ESU-BP/BEL (5.0 m)



K.Storz MIS Instrument Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BF (3.0 m), ESU-BP/BFL (5.0 m)



Ethicon/Everest Forceps Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BG (3.0 m), ESU-BP/BGL (5.0 m)



Ethicon Power Star Scissors Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BH (3.0 m), ESU-BP/BHL (5.0 m)



Ethicon Power Star Scissors Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BHT (3.0 m), ESU-BP/BHTL (5.0 m)



Aesculap Bipojet Scissors Connector

Integra, Martin, Aesculap, Berchtold, Sutter
BM780, Wolf

REF: ESU-BP/BI (3.0 m), ESU-BP/BIL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

Instrument Side Connector ▼



● **BiTech Scissors Connector**

Integra, Martin, Aesculap, Berchtold, Sutter BM780, Wolf

REF: ESU-BP/BK (3.0 m), ESU-BP/BKL (5.0 m)



● **To-Bite & Calvian Forceps v& Coagulation Clamps Connector**

Integra, Martin, Aesculap, Berchtold, Sutter BM780, Wolf

REF: ESU-BP/BL (3.0 m), ESU-BP/BLL (5.0 m)



● **K.Storz Resectoscope Connector**

Integra, Martin, Aesculap, Berchtold, Sutter BM780, Wolf

REF: ESU-BP/BM (3.0 m), ESU-BP/BML (5.0 m)



● **Olympus HiQ+Hand Instrument Connector**

Integra, Martin, Aesculap, Berchtold, Sutter BM780, Wolf

A60002C

REF: ESU-BP/BN (3.0 m), ESU-BP/BNL (5.0 m)



● **Tontarra, Wolf Resectoscope Connector**

Integra, Martin, Aesculap, Berchtold, Sutter BM780, Wolf

REF: ESU-BP/BO(3.0 m), ESU-BP/BOL(5.0 m)



● **Euro Forceps Connector**

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CA (3.0 m), ESU-BP/CAL (5.0 m)



● **U.S. Forceps Connector (Angled)**

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CB (3.0 m), ESU-BP/CBL (5.0 m)



● **U.S. Forceps Connector**

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CC (3.0 m), ESU-BP/CCL (5.0 m)



● **Aesculap Forceps Connector**

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CD (3.0 m), ESU-BP/CDL (5.0 m)



● **R.Wolf M.I.S. Instrument Connector**

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CE (3.0 m), ESU-BP/CEL (5.0 m)



● **K.Storz M.I.S. Instrument Connector**

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CF (3.0 m), ESU-BP/CFL (5.0 m)

** Depends on Generator socket

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

▼ Instrument Side Connector



● Ethicon / Everest Forceps Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CG (3.0 m), ESU-BP/CGL (5.0 m)



● Ethicon Power Stars Scissors Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CH (3.0 m), ESU-BP/CHL (5.0 m)



● Ethicon Power Stars Scissors Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CHT (3.0 m), ESU-BP/CHTL (5.0 m)



● Aesculap Bipojet Scissors Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CI (3.0 m), ESU-BP/CIL (5.0 m)



● BiTech Scissors Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CK (3.0 m), ESU-BP/CKL (5.0 m)



● To-Bite & Calvian Forceps & Coagulation Clamps Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CL (3.0 m), ESU-BP/CLL (5.0 m)



● K.Storz Resectoscope Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50

REF: ESU-BP/CM (3.0 m), ESU-BP/CML (5.0 m)



● Olympus HiQ+Hand Instrument Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50, Olympus HF 300


A6000C
REF: ESU-BP/CN (3.0 m), ESU-BP/CNL (5.0 m)



● Tontarra, Wolf Resectoscope Connector

Erbe ACC/ICC, Erbe T Series, Erbe VIO **, Siemens, Sutter BM50


REF: ESU-BP/CO (3.0 m), ESU-BP/COL (5.0 m)



● Euro Forceps Connector

Erbe, EMC/Dolley

REF: ESU-BP/DA (3.0 m), ESU-BP/DAL (5.0 m)



● U.S. Forceps Connector (Angled)

Erbe, EMC/Dolley

REF: ESU-BP/DB (3.0 m), ESU-BP/DBL (5.0 m)

** Depends on Generator socket

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

▼ Instrument Side Connector

			<p>● U.S. Forceps Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DC (3.0 m), ESU-BP/DCL (5.0 m)</p>
			<p>● Aesculap Forceps Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DD (3.0 m), ESU-BP/DDDL (5.0 m)</p>
			<p>● R.Wolf M.I.S. Instrument Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DE (3.0 m), ESU-BP/DEL (5.0 m)</p>
			<p>● K.Storz M.I.S. Instrument Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DF (3.0 m), ESU-BP/DFL (5.0 m)</p>
			<p>● Ethicon / Everest Forceps Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DG (3.0 m), ESU-BP/DGL (5.0 m)</p>
			<p>● Ethicon Power Star Scissors Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DH (3.0 m), ESU-BP/DHL (5.0 m)</p>
			<p>● Ethicon Power Star Scissors Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DHT (3.0 m), ESU-BP/DHTL (5.0 m)</p>
			<p>● Aesculap Bipojet Scissors Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DI (3.0 m), ESU-BP/DIL (5.0 m)</p>
			<p>● BiTech Scissors Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DK (3.0 m), ESU-BP/DKL (5.0 m)</p>
			<p>● To-Bite & Calvian Forceps & Coagulation Clamps Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DL (3.0 m), ESU-BP/DLL (5.0 m)</p>

** Depends on Generator socket

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

▼ Instrument Side Connector

			<p>● K.Storz Resectoscope Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DM (3.0 m), ESU-BP/DML (5.0 m)</p>
			<p>● Olympus HiQ+Hand Instrument Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DN (3.0 m), ESU-BP/DNL (5.0 m)</p>
			<p>● Tontarra, Wolf Resectoscope Connector</p> <p>Erbe, EMC/Dolley</p> <p>REF: ESU-BP/DO (3.0 m), ESU-BP/DOL (5.0 m)</p>
			<p>● Euro Forceps Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FA (3.0 m), ESU-BP/FAL (5.0 m)</p>
			<p>● U.S. Forceps Connector (Angled)</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FB (3.0 m), ESU-BP/FBL (5.0 m)</p>
			<p>● U.S. Forceps Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FC (3.0 m), ESU-BP/FCL (5.0 m)</p>
			<p>● Aesculap Forceps Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FD (3.0 m), ESU-BP/FDL (5.0 m)</p>
			<p>● R.Wolf M.I.S. Instrument Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FE (3.0 m), ESU-BP/FEL (5.0 m)</p>
			<p>● K.Storz M.I.S. Instrument Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FF (3.0 m), ESU-BP/FFL (5.0 m)</p>
			<p>● Ethicon / Everest Forceps Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FG (3.0 m), ESU-BP/FGL (5.0 m)</p>

** Depends on Generator socket

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

▼ Instrument Side Connector

			<p>● Ethicon Power Star Scissors Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FH (3.0 m), ESU-BP/FHL (5.0 m)</p>
			<p>● Ethicon / Everest Forceps Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FHT (3.0 m), ESU-BP/FHTL (5.0 m)</p>
			<p>● Aesculap Bipojet Scissors Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FI (3.0 m), ESU-BP/FIL (5.0 m)</p>
			<p>● BiTech Scissors Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FK (3.0 m), ESU-BP/FKL (5.0 m)</p>
			<p>● To-Bite & Calvian Forceps & Coagulation Clamps Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FL (3.0 m), ESU-BP/FLL (5.0 m)</p>
			<p>● K.Storz Resectoscope Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FM (3.0 m), ESU-BP/FML (5.0 m)</p>
			<p>● Olympus HiQ+Hand Instrument Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus UES-20/30/40 A60003C</p> <p>REF: ESU-BP/FN (3.0 m), ESU-BP/FNL (5.0 m)</p>
			<p>● Tontarra, Wolf Resectoscope Connector</p> <p>Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus</p> <p>REF: ESU-BP/FO(3.0 m), ESU-BP/FOL(5.0 m)</p>
			<p>● Euro Forceps Connector</p> <p>Valleylab Autobipolar</p> <p>REF: ESU-BP/GA (3.0 m), ESU-BP/GAL (5.0 m)</p>


** Depends on Generator socket

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector


▼ Instrument Side Connector



● U.S. Forceps Connector

Valleylab Autobipolar


REF: ESU-BP/GC (3.0 m), ESU-BP/GCL (5.0 m)



● Aesculap Forceps Connector

Valleylab Autobipolar

REF: ESU-BP/GD (3.0 m), ESU-BP/GDL (5.0 m)



● R.Wolf M.I.S. Instrument Connector

Valleylab Autobipolar


REF: ESU-BP/GE (3.0 m), ESU-BP/GEL (5.0 m)



● K.Storz M.I.S. Instrument Connector

Valleylab Autobipolar


REF: ESU-BP/GF (3.0 m), ESU-BP/GFL (5.0 m)



● Ethicon / Everest Forceps Connector

Valleylab Autobipolar


REF: ESU-BP/GG (3.0 m), ESU-BP/GGL (5.0 m)



● Ethicon Power Star Scissors Connector

Valleylab Autobipolar


REF: ESU-BP/GH (3.0 m), ESU-BP/GHL (5.0 m)



● Ethicon Power Star Scissors Connector

Valleylab Autobipolar


REF: ESU-BP/GHT (3.0 m), ESU-BP/GHTL (5.0 m)



● Aesculap Bipojet Scissors Connector

Valleylab Autobipolar



REF: ESU-BP/GI (3.0 m), ESU-BP/GIL (5.0 m)



● BiTech Scissors Connector

Valleylab Autobipolar

REF: ESU-BP/GK (3.0 m), ESU-BP/GKL (5.0 m)



● To-Bite & Calvian Forceps & Coagulation Clamps Connector

Valleylab Autobipolar


REF: ESU-BP/GL (3.0 m), ESU-BP/GLL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector


▼ Instrument Side Connector



● **K.Storz Resectoscope Connector**

Valleylab Autobipolar


REF: ESU-BP/GM (3.0 m), ESU-BP/GML (5.0 m)



● **Euro Forceps Connector**

Valleylab Force Triad


REF: ESU-BP/RA (3.0 m), ESU-BP/RAL (5.0 m)



● **U.S. Forceps Connector (Angled)**

Valleylab Force Triad

REF: ESU-BP/RB (3.0 m), ESU-BP/RBL (5.0 m)



● **U.S. Forceps Connector**

Valleylab Force Triad


REF: ESU-BP/RC (3.0 m), ESU-BP/RCL (5.0 m)



● **Aesculap Forceps Connector**

Valleylab Force Triad


REF: ESU-BP/RD (3.0 m), ESU-BP/RDL (5.0 m)



● **R.Wolf M.I.S. Instrument Connector**

Valleylab Force Triad


REF: ESU-BP/RE (3.0 m), ESU-BP/REL (5.0 m)



● **K.Storz M.I.S. Instrument Connector**

Valleylab Force Triad


REF: ESU-BP/RF (3.0 m), ESU-BP/RFL (5.0 m)



● **Ethicon / Everest Forceps Connector**

Valleylab Force Triad


REF: ESU-BP/RG (3.0 m), ESU-BP/RGL (5.0 m)



● **Ethicon Power Star Scissors Connector**

Valleylab Force Triad

REF: ESU-BP/RH (3.0 m), ESU-BP/RHL (5.0 m)



● **Ethicon Power Star Scissors Connector**

Valleylab Force Triad































REF: ESU-BP/RHT (3.0 m), ESU-BP/RHTL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

▼ Instrument Side Connector




























			<p>● Aesculap Bipojet Scissors Connector Valleylab Force Triad REF: ESU-BP/RI (3.0 m), ESU-BP/RIL (5.0 m)</p>
			<p>● BiTech Scissors Connector Valleylab Force Triad REF: ESU-BP/RK (3.0 m), ESU-BP/RKL (5.0 m)</p>
			<p>● To-Bite & Calvian Forceps & Coagulation Clamps Connector Valleylab Force Triad REF: ESU-BP/RL (3.0 m), ESU-BP/RLL (5.0 m)</p>
			<p>● K.Storz Resectoscope Connector Valleylab Autobipolar REF: ESU-BP/RM (3.0 m), ESU-BP/RML (5.0 m)</p>
			<p>● Euro Forceps Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HA (3.0 m), ESU-BP/HAL (5.0 m)</p>
			<p>● U.S. Forceps Connector (Angled) Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HB (3.0 m), ESU-BP/HBL (5.0 m)</p>
			<p>● U.S. Forceps Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HC (3.0 m), ESU-BP/HCL (5.0 m)</p>
			<p>● Aesculap Forceps Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HD (3.0 m), ESU-BP/HDL (5.0 m)</p>
			<p>● R.Wolf M.I.S. Instrument Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HE (3.0 m), ESU-BP/HEL (5.0 m)</p>
			<p>● K.Storz M.I.S. Instrument Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HF (3.0 m), ESU-BP/HFL (5.0 m)</p>

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

▼ Instrument Side Connector


 <p>2.4 mm</p>		 <p>2.0 mm</p> <p>● Ethicon / Everest Forceps Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HG (3.0 m), ESU-BP/HGL (5.0 m)</p>
 <p>2.4 mm</p>		 <p>2.5 mm</p> <p>● Ethicon Power Star Scissors Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HH (3.0 m), ESU-BP/HHL (5.0 m)</p>
 <p>2.4 mm</p>		 <p>2.0 mm</p> <p>● Aesculap Bipojet Scissors Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HI (3.0 m), ESU-BP/HIL (5.0 m)</p>
 <p>2.4 mm</p>		 <p>4.2 mm 8.0 mm</p> <p>● BiTech Scissors Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HK (3.0 m), ESU-BP/HKL (5.0 m)</p>
 <p>2.4 mm</p>		 <p>2 mm 4.4 mm</p> <p>● To-Bite & Calvian Forceps & Coagulation Clamps Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HL (3.0 m), ESU-BP/HLL (5.0 m)</p>
 <p>2.4 mm</p>		 <p>6.5 mm</p> <p>● K.Storz Resectoscope Connector Aesculap, GK4, GK20, GK21, GK30, GK34, GK50, GK51 REF: ESU-BP/HM (3.0 m), ESU-BP/HML (5.0 m)</p>
 <p>6.0 mm</p>		 <p>7.5 mm 4.0 mm</p> <p>● Euro Forceps Connector Alsa Excell MCDS Series, Alsatom SU/MPC Series REF: ESU-BP/JA (3.0 m), ESU-BP/JAL (5.0 m)</p>
 <p>6.0 mm</p>		 <p>1.8 mm</p> <p>● U.S. Forceps Connector (Angled) Alsa Excell MCDS Series, Alsatom SU/MPC Series REF: ESU-BP/JB (3.0 m), ESU-BP/JBL (5.0 m)</p>
 <p>6.0 mm</p>		 <p>1.8 mm</p> <p>● U.S. Forceps Connector Alsa Excell MCDS Series, Alsatom SU/MPC Series REF: ESU-BP/JC (3.0 m), ESU-BP/JCL (5.0 m)</p>

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector


▼ Instrument Side Connector



● **Aesculap Forceps Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JD (3.0 m), ESU-BP/JDL (5.0 m)



● **R.Wolf M.I.S. Instrument Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JE (3.0 m), ESU-BP/JEL (5.0 m)



● **K.Storz M.I.S. Instrument Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JF (3.0 m), ESU-BP/JFL (5.0 m)



● **Ethicon / Everest Forceps Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JG (3.0 m), ESU-BP/JGL (5.0 m)



● **Ethicon Power Star Scissors Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JH (3.0 m), ESU-BP/JHL (5.0 m)



● **Aesculap Bipojet Scissors Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JI (3.0 m), ESU-BP/JIL (5.0 m)



● **BiTech Scissors Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JK (3.0 m), ESU-BP/JKL (5.0 m)



● **To-Bite & Calvin Forceps v& Coagulation Clamps Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JL (3.0 m), ESU-BP/JLL (5.0 m)



● **K.Storz Resectoscope Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JM (3.0 m), ESU-BP/JML (5.0 m)



● **Tontarra, Wolf Resectoscope Connector**

Also Excell MCDS Series, Alsatom SU/MPC Series

REF: ESU-BP/JO (3.0 m), ESU-BP/JOL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

▼ Instrument Side Connector



● Euro Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KA (3.0 m), ESU-BP/KAL (5.0 m)



● Euro Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KAT (3.0 m), ESU-BP/KATL (5.0 m)



● U.S. Forceps Connector (Angled)

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KB (3.0 m), ESU-BP/KBL (5.0 m)



● U.S. Forceps Connector (Angled)

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KBT (3.0 m), ESU-BP/KBTL (5.0 m)



● U.S. Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KC (3.0 m), ESU-BP/KCL (5.0 m)



● U.S. Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KCT (3.0 m), ESU-BP/KCTL (5.0 m)



● Aesculap Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KD (3.0 m), ESU-BP/KDL (5.0 m)



● Aesculap Forceps Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KDT (3.0 m), ESU-BP/KDTL (5.0 m)



● R.Wolf MIS Instrument Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KE (3.0 m), ESU-BP/KEL (5.0 m)



● R.Wolf MIS Instrument Connector

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KET (3.0 m), ESU-BP/KETL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

▼ Generator Side Connector

Instrument Side Connector ▼



● **K.Storz M.I.S. Instrument Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab


REF: ESU-BP/KF (3.0 m), ESU-BP/KFL (5.0 m)



● **K.Storz M.I.S. Instrument Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KFT (3.0 m), ESU-BP/KFTL (5.0 m)



● **Ethicon/Everest Forceps Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KG (3.0 m), ESU-BP/KGL (5.0 m)



● **Ethicon/Everest Forceps Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KGT (3.0 m), ESU-BP/KGTL (5.0 m)



● **Ethicon Power Star Scissors Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab


REF: ESU-BP/KH (3.0 m), ESU-BP/KHL (5.0 m)



● **Aesculap Bipojet Scissors Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab


REF: ESU-BP/KI (3.0 m), ESU-BP/KIL (5.0 m)



● **Artery Sealer Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab


REF: ESU-BP/KJ (3.0 m), ESU-BP/KJL (5.0 m)



● **BiTech Scissors Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KK (3.0 m), ESU-BP/KKL (5.0 m)



● **To-Bite & Calvian Forceps & Coagulation Clamps Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KL (3.0 m), ESU-BP/KLL (5.0 m)



● **K.Storz Resectoscope Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KM(3.0 m), ESU-BP/KML (5.0 m)



● **Tontarra, Wolf Resectoscope Connector**

U.S. Standard, Erbe ICC, Alsa MCD Series, Eschmann, Petaş, Conmed, Valleylab

REF: ESU-BP/KO(3.0 m), ESU-BP/KOL(5.0 m)

Electrosurgery Cable & Accessories

High Frequency Bipolar Cables

4 mm adjustable connector with double banana for K series bipolar cables



Eschmann
18 mm



Erbe
22 mm



Valleylab
28 mm



ConMed
30 mm

▼ Generator Side Connector

▼ Instrument Side Connector



● **Olympus Working Element A220 60B/61B Connector**
Olympus UES-40
WA00013A
REF: ESU-BP/OWE(4.0 m)



● **Tontarra, Wolf Resectoscope Connector**
Olympus UES-40

REF: ESU-BP/OO(3.0 m), ESU-BP/OOL(5.0 m)



● **Olympus Working Element A220 60B/61B Connector**

Bowa, EMC, Lamidey, Erbe VIO **, Valleylab, Ethicon, Olympus


REF: ESU-BP/FWE(4.0 m)



● **Olympus Working Element A220 60B/61B Connector**

Integra, Martin, Aesculap, Berchtold, Sutter BM780, Wolf


REF: ESU-BP/BWE(4.0 m)



● **Tontarra, Wolf Resectoscope Connector**

Olympus / Gyrus PK-SP Generator

REF: ESU-BP/LO(4.0 m)



● **Olympus Working Element A220 60B/61B Connector**

Olympus / Gyrus PK-SP Generator

REF: ESU-BP/LWE(4.0 m)



● **Conmed Handswitching Forceps Connector**

Conmed System 2450, System 5000

REF: ESU-BP/CONA (3.0 m), ESU-BP/CONAL (5.0 m)



● **Gyrus ACMI Surgical EIWE-BRPK Iglesias PK Working Element Connector**

Gyrus ACMI PlasmaKinetic® SuperPulse Generator


REF: ESU-BP/GYR (3.0 m), ESU-BP/GYRL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Adapters



▼ Instrument Side Connector

▼ Generator Side Connector




Connecting finger switch handles with 3 pin plugs, for **Erbe T Series** generators.

REF: ESU-ADP/01



Connecting finger switch handles with 3 pin plugs, for **Erbe ACC / ICC, Erbe VIO** generators.

REF: ESU-ADP/02



Connecting finger switch handles with 3 pin plugs for **Martin, Integra** generators.

REF: ESU-ADP/03




4 mm banana to 6 mm Alsace Excel connector adapter for **Alsace Excel** generator

REF: ESU-ADP/04 (0.3 m)



4 mm banana to Martin connector adapter for **KLS Martin** generator

REF: ESU-ADP/05 (0.25 m)



Electrosurgery Cable & Accessories

High Frequency Monopolar Cables

Instrument Side Connector

Generator Side Connector



• Erbe T Series, Martin, Integra, Petaş
 Ø4.0 mm Banana Laparoscope / Endoscope cable

REF: ESU-MP/AE (3.0 m), ESU-MP/AEL (5.0 m)



• U.S. Standard, Erbe T Series, EschmannTD300 - TD311 - DS302, Sutter BM780, Martin, Integra, Aesculap, Berchtold
 Ø4.0 mm Female Laparoscope / Endoscope cable

REF: ESU-MP/BE (3.0 m), ESU-MP/BEL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab
 Ø4.0 mm Female Laparoscope / Endoscope cable

REF: ESU-MP/BF (3.0 m), ESU-MP/BFL (5.0 m)



• Aesculap, Berchtold, Erbe T Series, Huttinger, Martin, Micromed, Siemens
 Ø4.0 mm Female Laparoscope / Endoscope cable

REF: ESU-MP/BG (3.0 m), ESU-MP/BGL (5.0 m)



• ACC Standard, ICC, Erbe VIO, Karl Storz
 Ø4.0 mm Female Laparoscope / Endoscope cable

REF: ESU-MP/BD (3.0 m), ESU-MP/BDL (5.0 m)



• Erbe / Erbe T Series
 Ø4.0 mm Female Laparoscope / Endoscope cable

REF: ESU-MP/BR (3.0 m), ESU-MP/BRL (5.0 m)





• U.S. Standard, Erbe T Series, Eschmann TD300,TD311,DS302, Sutter BM780, Martin, Integra, Aesculap, Berchtold
 Ø4.0 mm Protected Banana Laparoscope / Endoscope cable

REF: ESU-MP/CE (3.0 m), ESU-MP/CEL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab
 Ø4.0 mm Protected Banana Laparoscope / Endoscope cable

REF: ESU-MP/CF (3.0 m), ESU-MP/CFL (5.0 m)



• Aesculap, Berchtold, Erbe T Series, Huttinger, Martin, Micromed, Siemens
 Ø4.0 mm Protected Banana Laparoscope / Endoscope cable

REF: ESU-MP/CG (3.0 m), ESU-MP/CGL (5.0 m)



• ACC Standard, ICC, Erbe VIO, Karl Storz
 Ø4.0 mm Protected Banana Laparoscope / Endoscope cable

REF: ESU-MP/CD (3.0 m), ESU-MP/CDL (5.0 m)



• Erbe / Erbe T Series
 Ø4.0 mm Protected Banana Laparoscope / Endoscope cable

REF: ESU-MP/CR (3.0 m), ESU-MP/CRL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Monopolar Cables

▼ Instrument Side Connector

Generator Side Connector ▼



• Erbe T Series, Eschmann TD300 TD311 DS302, Sutter BM780, Martin, Integra, Aesculap, Berchtold

Ø4 mm female laparoscope/endoscope cable

REF: ESU-MP/ DE (3.0 m), ESU-MP/ DEL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab

Ø4 mm female laparoscope/endoscope cable

REF: ESU-MP/ DF (3.0 m), ESU-MP/ DFL (5.0 m)



• Aesculap, Berchtold, Erbe T Series, Huttinger, Martin, Micromed, Siemens

Ø4 mm female laparoscope/endoscope cable

REF: ESU-MP/ DG (3.0 m), ESU-MP/ DGL (5.0 m)



• ACC Standard, ICC, Erbe VIO, Karl Storz

Ø4 mm female laparoscope/endoscope cable

REF: ESU-MP/ DD (3.0 m), ESU-MP/ DDL (5.0 m)



• U.S. Standard, Erbe T Series, Eschmann TD300 TD311 DS302, Sutter BM780, Martin, Integra, Aesculap, Berchtold

Ø2.0 mm Comeg / Wolf / Storz Resectoscope cable

REF: ESU-MP/ JE (3.0 m), ESU-MP/ JEL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab

Ø2.0 mm Comeg / Wolf / Storz Resectoscope cable

REF: ESU-MP/ JF (3.0 m), ESU-MP/ JFL (5.0 m)



• Aesculap, Berchtold, Erbe T Series, Huttinger, Martin, Micromed, Siemens

Ø2.0 mm Comeg / Wolf / Storz Resectoscope cable

REF: ESU-MP/ JG (3.0 m), ESU-MP/ JGL (5.0 m)



• ACC Standard, ICC, Erbe VIO, Karl Storz

Ø2.0 mm Comeg / Wolf / Storz Resectoscope cable

REF: ESU-MP/ JD (3.0 m), ESU-MP/ JDL (5.0 m)



• Erbe / Erbe T Series

Ø2.0 mm Comeg / Wolf / Storz Resectoscope cable

REF: ESU-MP/ JR (3.0 m), ESU-MP/ JRL (5.0 m)



• Erbe T Series, Eschmann TD300 TD311 DS302, Sutter BM780, Martin, Integra, Aesculap, Berchtold

Storz Resectoscope cable




























REF: ESU-MP/ KE (3.0 m), ESU-MP/ KEL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Monopolar Cables

▼ Instrument Side Connector

▼ Generator Side Connector

			<ul style="list-style-type: none"> • Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab Storz Resectoscope cable <p>REF: ESU-MP/KF (3.0 m), ESU-MP/KFL (5.0 m)</p>
			<ul style="list-style-type: none"> • Aesculap, Berchtold, Erbe T Series, Huttering, Martin, Micromed, Siemens Storz Resectoscope cable <p>REF: ESU-MP/KG (3.0 m), ESU-MP/KGL (5.0 m)</p>
			<ul style="list-style-type: none"> • ACC Standard, ICC, Erbe VIO, Karl Storz Storz Resectoscope cable <p>REF: ESU-MP/KD (3.0 m), ESU-MP/KDL (5.0 m)</p>
			<ul style="list-style-type: none"> • Erbe / Erbe T Series Storz Resectoscope cable <p>REF: ESU-MP/KR (3.0 m), ESU-MP/KRL (5.0 m)</p>
			<ul style="list-style-type: none"> • Erbe T Series, Eschmann TD300, TD311, DS302, Sutter BM780, Martin, Integra, Aesculap, Berchtold Ø3.0 mm Female Olympus - Fujinon Endoscope Instrument cable <p>REF: ESU-MP/LE (3.0 m), ESU-MP/LEL (5.0 m)</p>
			<ul style="list-style-type: none"> • Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab Ø3.0 mm Female Olympus - Fujinon Endoscope Instrument cable <p>REF: ESU-MP/LF (3.0 m), ESU-MP/LFL (5.0 m)</p>
			<ul style="list-style-type: none"> • Aesculap, Berchtold, Erbe T Series, Huttering, Martin, Micromed, Siemens Ø3.0 mm Female Olympus - Fujinon Endoscope Instrument cable <p>REF: ESU-MP/LG (3.0 m), ESU-MP/LGL (5.0 m)</p>
			<ul style="list-style-type: none"> • ACC Standard, ICC, Erbe VIO, Karl Storz Ø3.0 mm Female Olympus - Fujinon Endoscope Instrument cable <p>REF: ESU-MP/LD (3.0 m), ESU-MP/LDL (5.0 m)</p>
			<ul style="list-style-type: none"> • Erbe ACC/ Erbe T Series Ø3.0 mm Female Olympus - Fujinon Endoscope Instrument cable <p>REF: ESU-MP/LR (3.0 m), ESU-MP/LRL (5.0 m)</p>

Electrosurgery Cable & Accessories

High Frequency Monopolar Cables

▼ Instrument Side Connector

▼ Generator Side Connector



• Olympus
 Ø3.0 mm Female Olympus - Fujinon Endoscope Instrument cable

REF: ESU-MP/LO (3.0 m), ESU-MP/LOL (5.0 m)



• Erbe T Series, Eschmann TD300 - TD311 - DS302 - Sutter BM780, Martin, Integra, Aesculap, Berchtold
 Ø4.8 mm UK Type Monopolar Forceps cable

REF: ESU-MP/ME (3.0 m), ESU-MP/MEL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab
 Ø4.8 mm UK Type Monopolar Forceps cable

REF: ESU-MP/MF (3.0 m), ESU-MP/MFL (5.0 m)



• Aesculap, Berchtold, Erbe T Series, Huttinger, Martin, Micromed, Siemens
 Ø4.8 mm UK Type Monopolar Forceps cable

REF: ESU-MP/MG (3.0 m), ESU-MP/MGL (5.0 m)



• ACC Standard, ICC, Erbe VIO, Karl Storz
 Ø4.8 mm UK Type Monopolar Forceps cable

REF: ESU-MP/MD (3.0 m), ESU-MP/MDL (5.0 m)



• Erbe ACC/ Erbe T Series
 Ø4.8 mm UK Type Monopolar Forceps cable

REF: ESU-MP/MR (3.0 m), ESU-MP/MRL (5.0 m)



• Erbe T Series, Eschmann TD300 - TD311 - DS302 - Sutter BM780, Martin, Integra, Aesculap, Berchtold
 Ø4.0 mm Socket Hexagonal Storz Instrument cable

REF: ESU-MP/NE (3.0 m), ESU-MP/NEL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab
 Ø4.0 mm Socket Hexagonal Storz Instrument cable

REF: ESU-MP/NF (3.0 m), ESU-MP/NFL (5.0 m)



• Aesculap, Berchtold, Erbe T Series, Huttinger, Martin, Micromed, Siemens
 Ø4.0 mm Socket Hexagonal Storz Instrument cable

REF: ESU-MP/NG (3.0 m), ESU-MP/NGL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Monopolar Cables

▼ Instrument Side Connector

▼ Generator Side Connector



• ACC Standard, ICC, Erbe VIO, Karl Storz
 Ø4.0 mm Socket Hexagonal Storz Instrument cable

REF: ESU-MP/ ND (3.0 m), ESU-MP/ NDL (5.0 m)



• Erbe ACC/ Erbe T Series
 Ø4.0 mm Socket Hexagonal Storz Instrument cable

REF: ESU-MP/ NR (3.0 m), ESU-MP/ NRL (5.0 m)



• Erbe T Series, Eschmann TD300 TD311 DS302, Sutter BM780, Martin, Integra, Aesculap, Berchtold
 Storz Pediatric Resectoscope cable

REF: ESU-MP/ PE (3.0 m), ESU-MP/ PEL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab
 Storz Pediatric Resectoscope cable

REF: ESU-MP/ PF (3.0 m), ESU-MP/ PFL (5.0 m)



• Aesculap, Berchtold, Erbe T Series, Huttering, Martin, Micromed, Siemens
 Storz Pediatric Resectoscope cable

REF: ESU-MP/ PG (3.0 m), ESU-MP/ PGL (5.0 m)



• ACC Standard, ICC, Erbe VIO, Karl Storz
 Storz Pediatric Resectoscope cable

REF: ESU-MP/ PD (3.0 m), ESU-MP/ PDL (5.0 m)



• Erbe ACC/ Erbe T Series
 Storz Pediatric Resectoscope cable

REF: ESU-MP/ PR (3.0 m), ESU-MP/ PRL (5.0 m)



• Erbe T Series, Eschmann TD300 - TD311 - DS302 - Sutter BM780, Martin, Integra, Aesculap, Berchtold
 Stryker Working Element Cable

REF: ESU-MP/ RE (3.0 m), ESU-MP/ REL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab
 Stryker Working Element Cable

REF: ESU-MP/ RF (3.0 m), ESU-MP/ RFL (5.0 m)

Electrosurgery Cable & Accessories

High Frequency Monopolar Cables

▼ Instrument Side Connector

▼ Generator Side Connector



• Aesculap, Berchtold, Erbe T Series, Huttlinger, Martin, Micromed, Siemens
Stryker Working Element Cable

REF: ESU-MP/ RG (3.0 m), ESU-MP/ RGL (5.0 m)



• ACC Standard, ICC, Erbe VIO, Karl Storz
Stryker Working Element Cable

REF: ESU-MP/ RD (3.0 m), ESU-MP/ RDL (5.0 m)



• Erbe ACC/ Erbe T Series
Stryker Working Element Cable

REF: ESU-MP/ RR (3.0 m), ESU-MP/ RRL (5.0 m)



• Erbe T Series, Eschmann TD300 - TD311 - DS302 - Sutter BM780, Martin, Integra, Aesculap, Berchtold
Olympus Resectoscope Cable A0392

REF: ESU-MP/ SE (3.0 m), ESU-MP/ SEL (5.0 m)



• Bard, Birtcher, Bovie, Bowa, Conmed, Eschmann, Medtrex, Valleylab, Olympus UES-30/40
Olympus Resectoscope Cable A0393

REF: ESU-MP/ SF (3.0 m), ESU-MP/ SFL (5.0 m)



• Aesculap, Berchtold, Erbe T Series, Huttlinger, Martin, Micromed, Siemens
Olympus Resectoscope Cable

REF: ESU-MP/ SG (3.0 m), ESU-MP/ SGL (5.0 m)



• ACC Standard, ICC, Erbe VIO, Karl Storz
Olympus Resectoscope Cable

REF: ESU-MP/ SD (3.0 m), ESU-MP/ SDL (5.0 m)



• Erbe ACC/ Erbe T Series
Olympus Resectoscope Cable

REF: ESU-MP/ SR (3.0 m), ESU-MP/ SRL (5.0 m)



• Olympus UES-10/20, Valleylab (Old)
Olympus Resectoscope cable A0391

REF: ESU-MP/ SO (3.0 m), ESU-MP/ SOL (5.0 m)



Reusable Dispersive Cables, Blue Series (Non-Autoclavable)

▼ Dispersive Pad Clamp (23.5mm)

▼ Generator Side Connector



Alsa, Bard, Berchtold, Bovie-Aaron 950 / 1200, Eschmann TD400/ TD800, Martin, Meditom

REF: RDC-A3 (3.0 m), RDC-A5 (5.0 m)



Erbe ACC/ICC, Erbe T Series, Erbe VIO, Martin

REF: RDC-B3 (3.0 m), RDC-B5 (5.0 m)



Valleylab REM, Bovie-Aaron 1250 / 2250 / 3250, Aspen, Bowa, Conmed, Petaş 500 / 600.

REF: RDC-C3 (3.0 m), RDC-C5 (5.0 m)



Valleylab NON-REM, Bowa, Conmed

REF: RDC-D3 (3.0 m), RDC-D5 (5.0 m)



Valleylab (Old Generation)

REF: RDC-E3 (3.0 m), RDC-E5 (5.0 m)



Bovie

REF: RDC-F3 (3.0 m), RDC-F5 (5.0 m)



Actoll, Diadon, Downs 450P2, Eschmann TD300 / TD311 / DS302

REF: RDC-G3 (3.0 m), RDC-G5 (5.0 m)



Petaş 400/ 450 (Old Generation), Eschmann (Adjustable)

REF: RDC-H3 (3.0 m), RDC-H5 (5.0 m)



Reusable Dispersive Cables, Black Series (Autoclavable)

▼ Dispersive Pad Clamp (23.5 mm)

▼ Generator Side Connector



Alsa, Bard, Berchtold, Bovie-Aaron 950 / 1200, Eschmann TD400/ TD800, Martin, Meditom

REF: RDC-A3A (3.0 m), RDC-A5A (5.0 m)



Erbe ACC/ICC, Erbe T Series, Erbe VIO, Martin

REF: RDC-B3A (3.0 m), RDC-B5A (5.0 m)



Valleylab REM, Bovie-Aaron 1250 / 2250 / 3250, Aspen, Bowa, Conmed, Petaş 500 / 600.

REF: RDC-C3A (3.0 m), RDC-C5A (5.0 m)



Valleylab NON-REM, Bowa, Conmed

REF: RDC-D3A (3.0 m), RDC-D5A (5.0 m)



Valleylab (Old Generation)

REF: RDC-E3A (3.0 m), RDC-E5A (5.0 m)



Bovie

REF: RDC-F3A (3.0 m), RDC-F5A (5.0 m)



Actoll, Diadon, Downs 450P2, Eschmann TD300 / TD311 / DS302

REF: RDC-G3A (3.0 m), RDC-G5A (5.0 m)



Petaş 400/ 450 (Old Generation), Eschmann (Adjustable)

REF: RDC-H3A (3.0 m), RDC-H5A (5.0 m)

Reusable Valleylab Style Dispersive Cables, Blue Series (Non-Autoclavable)

▼ Dispersive Pad Clamp (32mm)

▼ Generator Side Connector



Valleylab REM, Bovie-Aaron 1250 / 2250 / 3250, Aspen, Bowa, Conmed, Petaş 500 / 600.

REF: RDC-VL3 (3.0 m), RDC-VL5 (5.0 m)



Erbe ACC/ICC, Erbe T Series, Erbe VIO, Martin

REF: RDC-VB3 (3.0 m), RDC-VB5 (5.0 m)



Alsa, Bard, Berchtold, Bovie-Aaron 950 / 1200, Eschmann TD400 / TD800, Martin, Meditom

REF: RDC-VA3 (3.0 m), RDC-VA5 (5.0 m)

Reusable Valleylab Style Dispersive Cables, Black Series (Autoclavable)

▼ Dispersive Pad Clamp (32mm)

▼ Generator Side Connector



Valleylab REM, Bovie-Aaron 1250 / 2250 / 3250, Aspen, Bowa, Conmed, Petaş 500 / 600

REF: RDC-VL3A (3.0 m), RDC-VL5A (5.0 m)



Erbe ACC/ICC, Erbe T Series, Erbe VIO, Martin

REF: RDC-VB3A (3.0 m), RDC-VB5A (5.0 m)



Alsa, Bard, Berchtold, Bovie-Aaron 950 / 1200, Eschmann TD400 / TD800, Martin, Meditom

REF: RDC-VA3A (3.0 m), RDC-VA5A (5.0 m)



Electrosurgery Cable & Accessories

Connection Cable for MEGADYNE Universal Patient Return Electrodes



● **FMT-MD1/24 (M2K-1)** Reusable Compatibility Cable, Standard Connector, 2.4m (8.0 ft) for Megadyne Mega Soft, Mega Soft Dual, Mega Soft Pediatric



● **FMT-MD2/44 (M2K-2)** Reusable Compatibility Cable, Standard Connector, 4.4m (14.4 ft) for Megadyne Mega Soft, Mega Soft Dual, Mega Soft Pediatric



● **FMT-MD5/44 (M2K-5)** Reusable Compatibility Cable, Extended 6.35mm Phone Plug, 4.4m (14.4 ft) for Megadyne Mega Soft, Mega Soft Dual, Mega Soft Pediatric



● **FMT-MD7 (M2K-7)** Replacement Pigtail Cable 0.82m (2.7 ft) for Megadyne Mega Soft, Mega Soft Dual, Mega Soft Pediatric



● **FMT-MD8/24 (M2K-8)** Reusable Compatibility Cable, Standard Connector, 2.4m (8.0 ft) for Megadyne Mega2000



● **FMT-MD9/44 (M2K-9)** Reusable Compatibility Cable, Standard Connector, 4.4m (14.4 ft) for Megadyne Mega2000



■ FMT foot control ESU pencil without cable

REF: FSP-24

Ø2.40 mm



Ø4.00 mm banana connector



REF: FSP-40


Ø4.00 mm



Ø4.00 mm banana connector



- Reusable Autoclavable (134°C / 18 min.), 100 uses minimum, non sterile.
- Polyacetal (POM) Material, 14.5 cm length.
- Machine Wash (operating room or laboratory).
- Compatible with all generator types.
- Accepts all standart 2.40 mm (3/32 inch) or 4.00 mm (5/32 inch) shaft electrodes.
- Without cable and steel blade.



FMT is registered trademark of **Metko Ltd.**



Headquarter & Factory

İvedik O.S.B. Ağaç İşleri Sanayi Sit.
1354. Cad. 1358.Sok. No:9
06378 Yenimahalle

ANKARA / TÜRKİYE

Tel (pbx) : +90 312 387 12 46
Fax : +90 312 387 12 51

Branch

Atatürk Mah. Sedef Cad. 38.Ada
Ata Plaza 3-5 Blokları Kat:2 No:57
34758 Ataşehir

İSTANBUL / TÜRKİYE

Tel (pbx) : +90 216 455 91 21
Fax : +90 216 455 48 24

www.metkold.com

Other products, logos and company names mentioned herein may be trademarks of their respective owners and are not to be taken as an endorsement of affiliation with FMT®. The information in this catalogue has been carefully checked and it is believed to be accurate. Images in this catalogue are representative of items offered. Actual product may have slight visual differences. In the interest of continued product development, Metko reserves the right to make changes and improvements to this catalogue and the product it described any time, without notice or obligation.